



VC Test

VC test flow



VC test comparison screen



VC test curve screen

VC Test Result			
Para.	Pred.	Meas.	%Pred.
VC-----[L]:	3.508	2.674	76.2%
TV-----[L]:		0.953	
ERV-----[L]:		0.339	
IRV-----[L]:		1.380	
IC-----[L]:		2.334	

[START]Retest, [←]Review, [↵]Enter

VC test results screen

VC test report

Patient Test Data		ID: 000000015 SVC		2009/08/27 14:04		
Name: _____		Para	Unit	Pred.	Meas.	%Pred.
2009/08/27 14:04		VC	L	4.613	3.767	81.6%
ID No.: 000000015		TV	L		0.501	
Gender: Male	Age: 30 yrs	ERV	L		1.494	
Height: 170 cm	Weight: 70 kg	IRV	L		1.772	
Race: Chinese		IC	L		2.273	
Temp.: 29 C	Humi.: 40 %					
Atoms: 760 mmHg						
Test Result:						


Doctor: _____						

[SVC]


VC test report

MVV Test


MVV test flow



MVV test comparison screen



MVV Test curve and results MVV test report




MVV test report

FVC test

FVC test flow

<div style="text-align: right; font-size: small;">2009/03/27 11:57:31</div> <p style="text-align: center;">Pulmonary Function Test</p> <ol style="list-style-type: none"> 1. User Info Entry 2. VC Testing 3. FVC Testing 4. MVV Testing 5. BronChodilator Test 6. UTILITY <div style="text-align: center; font-size: x-small;">[▲][▼]Select, [↵]Enter</div>	<div style="text-align: right; font-size: small;">2009/03/27 11:57:31</div> <p style="text-align: center;">FVC Test</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Index</th> <th style="text-align: left;">FVC (L)</th> <th style="text-align: left;">FEV1 (L)</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">1</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">2</td> <td></td> <td></td> </tr> <tr> <td style="text-align: center;">3</td> <td></td> <td></td> </tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[0]Save, [↵]Enter, [↶]Return, [START]</div>	Index	FVC (L)	FEV1 (L)	1			2			3		
Index	FVC (L)	FEV1 (L)											
1													
2													
3													

FVC test comparison screen



FVC test curve

<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="4" style="text-align: center;">FVC Test Result</th> <th style="text-align: right;">1/3</th> </tr> <tr> <th style="text-align: left;">Para.</th> <th style="text-align: left;">Pred.</th> <th style="text-align: left;">Meas.</th> <th style="text-align: left;">%Pred.</th> <th></th> </tr> </thead> <tbody> <tr> <td>FVC</td><td>[L]: 6.569</td><td>3.240</td><td>49.3%</td><td></td></tr> <tr> <td>FEV0.5</td><td>[L]: 3.917</td><td>2.788</td><td>71.1%</td><td></td></tr> <tr> <td>FEV1.0</td><td>[L]: 5.684</td><td>3.070</td><td>54.0%</td><td></td></tr> <tr> <td>FEV3.0</td><td>[L]: 6.348</td><td>3.070</td><td>48.4%</td><td></td></tr> <tr> <td>FEV1.0/G</td><td>[%]: 84.3%</td><td>94.7%</td><td>112.3%</td><td></td></tr> <tr> <td>FEV1.0/T</td><td>[%]:</td><td>0.0%</td><td></td><td></td></tr> <tr> <td>FEV3.0/G</td><td>[%]:</td><td>94.7%</td><td></td><td></td></tr> <tr> <td>FEV3.0/T</td><td>[%]:</td><td>0.0%</td><td></td><td></td></tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[↶]Return, [▼]Next</div>	FVC Test Result				1/3	Para.	Pred.	Meas.	%Pred.		FVC	[L]: 6.569	3.240	49.3%		FEV0.5	[L]: 3.917	2.788	71.1%		FEV1.0	[L]: 5.684	3.070	54.0%		FEV3.0	[L]: 6.348	3.070	48.4%		FEV1.0/G	[%]: 84.3%	94.7%	112.3%		FEV1.0/T	[%]:	0.0%			FEV3.0/G	[%]:	94.7%			FEV3.0/T	[%]:	0.0%			<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="4" style="text-align: center;">FVC Test Result</th> <th style="text-align: right;">2/3</th> </tr> <tr> <th style="text-align: left;">Para.</th> <th style="text-align: left;">Pred.</th> <th style="text-align: left;">Meas.</th> <th style="text-align: left;">%Pred.</th> <th></th> </tr> </thead> <tbody> <tr> <td>Vent</td><td>[L]:</td><td>0.050</td><td></td><td></td></tr> <tr> <td>EX Time</td><td>[s]:</td><td>6.0</td><td></td><td></td></tr> <tr> <td>MMF</td><td>[L/s]:</td><td>6.270</td><td>6.943</td><td>110.7%</td></tr> <tr> <td>PEF</td><td>[L/s]:</td><td>11.368</td><td>7.098</td><td>62.4%</td></tr> <tr> <td>MEF75</td><td>[L/s]:</td><td>3.481</td><td>5.231</td><td>150.2%</td></tr> <tr> <td>MEF50</td><td>[L/s]:</td><td>7.126</td><td>7.098</td><td>99.6%</td></tr> <tr> <td>MEF25</td><td>[L/s]:</td><td>10.334</td><td>5.702</td><td>55.1%</td></tr> <tr> <td>FIVC</td><td>[L]:</td><td>3.304</td><td></td><td></td></tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[▲][▼]Page Up and Down</div>	FVC Test Result				2/3	Para.	Pred.	Meas.	%Pred.		Vent	[L]:	0.050			EX Time	[s]:	6.0			MMF	[L/s]:	6.270	6.943	110.7%	PEF	[L/s]:	11.368	7.098	62.4%	MEF75	[L/s]:	3.481	5.231	150.2%	MEF50	[L/s]:	7.126	7.098	99.6%	MEF25	[L/s]:	10.334	5.702	55.1%	FIVC	[L]:	3.304		
FVC Test Result				1/3																																																																																																	
Para.	Pred.	Meas.	%Pred.																																																																																																		
FVC	[L]: 6.569	3.240	49.3%																																																																																																		
FEV0.5	[L]: 3.917	2.788	71.1%																																																																																																		
FEV1.0	[L]: 5.684	3.070	54.0%																																																																																																		
FEV3.0	[L]: 6.348	3.070	48.4%																																																																																																		
FEV1.0/G	[%]: 84.3%	94.7%	112.3%																																																																																																		
FEV1.0/T	[%]:	0.0%																																																																																																			
FEV3.0/G	[%]:	94.7%																																																																																																			
FEV3.0/T	[%]:	0.0%																																																																																																			
FVC Test Result				2/3																																																																																																	
Para.	Pred.	Meas.	%Pred.																																																																																																		
Vent	[L]:	0.050																																																																																																			
EX Time	[s]:	6.0																																																																																																			
MMF	[L/s]:	6.270	6.943	110.7%																																																																																																	
PEF	[L/s]:	11.368	7.098	62.4%																																																																																																	
MEF75	[L/s]:	3.481	5.231	150.2%																																																																																																	
MEF50	[L/s]:	7.126	7.098	99.6%																																																																																																	
MEF25	[L/s]:	10.334	5.702	55.1%																																																																																																	
FIVC	[L]:	3.304																																																																																																			
<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="4" style="text-align: center;">FVC Test Result</th> <th style="text-align: right;">3/3</th> </tr> <tr> <th style="text-align: left;">Para.</th> <th style="text-align: left;">Pred.</th> <th style="text-align: left;">Meas.</th> <th style="text-align: left;">%Pred.</th> <th></th> </tr> </thead> <tbody> <tr> <td>FIV0.5</td><td>[L]:</td><td>1.714</td><td></td><td></td></tr> <tr> <td>FIV1.0</td><td>[L]:</td><td>3.304</td><td></td><td></td></tr> <tr> <td>FIV1.0/FVC</td><td>[%]:</td><td>101.9%</td><td></td><td></td></tr> <tr> <td>FIV1.0/FIVC</td><td>[%]:</td><td>100.0%</td><td></td><td></td></tr> <tr> <td>PIF</td><td>[L/s]:</td><td>7.386</td><td></td><td></td></tr> <tr> <td>MIF50%</td><td>[L/s]:</td><td>7.260</td><td></td><td></td></tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[START]Retest, [▲]Page Up, [↵]Enter, [0]Save, [PRINT]</div>	FVC Test Result				3/3	Para.	Pred.	Meas.	%Pred.		FIV0.5	[L]:	1.714			FIV1.0	[L]:	3.304			FIV1.0/FVC	[%]:	101.9%			FIV1.0/FIVC	[%]:	100.0%			PIF	[L/s]:	7.386			MIF50%	[L/s]:	7.260			<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2" style="text-align: center;">FVC Test Result</th> <th style="text-align: right;">4/4</th> </tr> </thead> <tbody> <tr> <td colspan="2">Diagnosis:</td> <td></td> </tr> <tr> <td colspan="2" style="text-align: center;">Lung diseases and degree:</td> <td></td> </tr> <tr> <td colspan="2" style="text-align: center;">Mild obstruction</td> <td></td> </tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[START]Retest, [▲]Page Up, [↵]Enter</div>	FVC Test Result		4/4	Diagnosis:			Lung diseases and degree:			Mild obstruction																																																		
FVC Test Result				3/3																																																																																																	
Para.	Pred.	Meas.	%Pred.																																																																																																		
FIV0.5	[L]:	1.714																																																																																																			
FIV1.0	[L]:	3.304																																																																																																			
FIV1.0/FVC	[%]:	101.9%																																																																																																			
FIV1.0/FIVC	[%]:	100.0%																																																																																																			
PIF	[L/s]:	7.386																																																																																																			
MIF50%	[L/s]:	7.260																																																																																																			
FVC Test Result		4/4																																																																																																			
Diagnosis:																																																																																																					
Lung diseases and degree:																																																																																																					
Mild obstruction																																																																																																					

FVC test results

⚠ Attention: Because one chooses different diagnosis mode in “Mode Setting” screen, the content and style of FVC test results will vary in fourth page.

Shown as in figure 4-16, the fourth page is the diagnosis results in ATS or NIOSH mode, if the mode is changed to DIAGNOSIS, the results is like right-hand figure.

FVC Test Result		4/4
Diagnosis:		
FEV1C	REST *	NORM
70%		
	COMB	OBST
		80% VC%
[START]Retest, [▲]Page Up, [↵]Enter		

In the figure: NORM means ordinary


REST means restricted type harm

OBST means obstruction type harm;


COMB means composite type harm (restraint + obstruction).

ult, you can't print test report at last.

FVC test report



One style of FVC report



Second style of FVC report

<div style="text-align: right; font-size: small;">2009/03/27 11:57:31</div> <p style="text-align: center;">Pulmonary Function Test</p> <ol style="list-style-type: none"> 1. User Info Entry 2. VC Testing 3. FVC Testing 4. MVV Testing 5. BronChodilator Test 6. UTILITY <div style="text-align: center; font-size: x-small;">[▲][▼]Select, [↵]Enter</div>	<p>It is in post-measurement!</p>
--	-----------------------------------

Post mode

<div style="text-align: right; font-size: small;">2009/03/27 11:57:31</div> <p style="text-align: center;">VC Test--Post</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Para.</th> <th style="text-align: center;">Pre.</th> <th style="text-align: center;">Post.</th> <th style="text-align: center;">%Pre.</th> </tr> </thead> <tbody> <tr> <td>VC-----[L]:</td> <td style="text-align: center;">2.700</td> <td style="text-align: center;">2.501</td> <td style="text-align: center;">92.6%</td> </tr> <tr> <td>TV-----[L]:</td> <td style="text-align: center;">0.511</td> <td style="text-align: center;">0.449</td> <td style="text-align: center;">87.8%</td> </tr> <tr> <td>ERV-----[L]:</td> <td style="text-align: center;">1.022</td> <td style="text-align: center;">0.94</td> <td style="text-align: center;">91.9%</td> </tr> <tr> <td>IRV-----[L]:</td> <td style="text-align: center;">1.167</td> <td style="text-align: center;">1.112</td> <td style="text-align: center;">95.2%</td> </tr> <tr> <td>IC-----[L]:</td> <td style="text-align: center;">1.678</td> <td style="text-align: center;">1.561</td> <td style="text-align: center;">93.0%</td> </tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[START]Retest, [←]Review, [↵]Enter</div>	Para.	Pre.	Post.	%Pre.	VC-----[L]:	2.700	2.501	92.6%	TV-----[L]:	0.511	0.449	87.8%	ERV-----[L]:	1.022	0.94	91.9%	IRV-----[L]:	1.167	1.112	95.2%	IC-----[L]:	1.678	1.561	93.0%	
Para.	Pre.	Post.	%Pre.																						
VC-----[L]:	2.700	2.501	92.6%																						
TV-----[L]:	0.511	0.449	87.8%																						
ERV-----[L]:	1.022	0.94	91.9%																						
IRV-----[L]:	1.167	1.112	95.2%																						
IC-----[L]:	1.678	1.561	93.0%																						

SVC results in post-measurement

<div style="text-align: right; font-size: small;">2009/03/27 11:57:31</div> <p style="text-align: center;">MVV Test--Post</p> <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left;">Para.</th> <th style="text-align: center;">Pre.</th> <th style="text-align: center;">Post.</th> <th style="text-align: center;">%Pre.</th> </tr> </thead> <tbody> <tr> <td>MVV---[L/min]:</td> <td style="text-align: center;">69.521</td> <td style="text-align: center;">61.061</td> <td style="text-align: center;">87.8%</td> </tr> <tr> <td>RR---[/min]:</td> <td style="text-align: center;">48</td> <td style="text-align: center;">57</td> <td style="text-align: center;">118.7%</td> </tr> <tr> <td>TV-----[L]:</td> <td style="text-align: center;">1.448</td> <td style="text-align: center;">1.071</td> <td style="text-align: center;">73.9%</td> </tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[START]Retest, [←]Review, [↵]Enter</div>	Para.	Pre.	Post.	%Pre.	MVV---[L/min]:	69.521	61.061	87.8%	RR---[/min]:	48	57	118.7%	TV-----[L]:	1.448	1.071	73.9%	
Para.	Pre.	Post.	%Pre.														
MVV---[L/min]:	69.521	61.061	87.8%														
RR---[/min]:	48	57	118.7%														
TV-----[L]:	1.448	1.071	73.9%														

MVV results in post-measurement



<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Para.</th> <th colspan="2">FVC</th> <th>1/4</th> </tr> <tr> <th>Pre.</th> <th>Post.</th> <th>%Pre.</th> </tr> </thead> <tbody> <tr> <td>FVC-----[L]:</td> <td style="text-align: center;">2.671</td> <td style="text-align: center;">2.493</td> <td style="text-align: center;">93.3%</td> </tr> <tr> <td>FEV0.5---[L]:</td> <td style="text-align: center;">1.753</td> <td style="text-align: center;">1.751</td> <td style="text-align: center;">99.8%</td> </tr> <tr> <td>FEV1.0---[L]:</td> <td style="text-align: center;">2.276</td> <td style="text-align: center;">2.247</td> <td style="text-align: center;">98.7%</td> </tr> <tr> <td>FEV3.0---[L]:</td> <td style="text-align: center;">2.616</td> <td style="text-align: center;">2.463</td> <td style="text-align: center;">94.1%</td> </tr> <tr> <td>FEV1.0/G---[%]:</td> <td style="text-align: center;">85.2%</td> <td style="text-align: center;">99.8%</td> <td style="text-align: center;">105.7%</td> </tr> <tr> <td>FEV1.0/T---[%]:</td> <td style="text-align: center;">84.2%</td> <td style="text-align: center;">89.8%</td> <td style="text-align: center;">106.6%</td> </tr> <tr> <td>FEV3.0/G---[%]:</td> <td style="text-align: center;">97.9%</td> <td style="text-align: center;">98.7%</td> <td style="text-align: center;">100.8%</td> </tr> <tr> <td>FEV3.0/T---[%]:</td> <td style="text-align: center;">96.8%</td> <td style="text-align: center;">98.4%</td> <td style="text-align: center;">101.6%</td> </tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[←]Return, [▼]Next</div>	Para.	FVC		1/4	Pre.	Post.	%Pre.	FVC-----[L]:	2.671	2.493	93.3%	FEV0.5---[L]:	1.753	1.751	99.8%	FEV1.0---[L]:	2.276	2.247	98.7%	FEV3.0---[L]:	2.616	2.463	94.1%	FEV1.0/G---[%]:	85.2%	99.8%	105.7%	FEV1.0/T---[%]:	84.2%	89.8%	106.6%	FEV3.0/G---[%]:	97.9%	98.7%	100.8%	FEV3.0/T---[%]:	96.8%	98.4%	101.6%	<table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th rowspan="2">Para.</th> <th colspan="2">FVC</th> <th>2/4</th> </tr> <tr> <th>Pre.</th> <th>Post.</th> <th>%Pre.</th> </tr> </thead> <tbody> <tr> <td>Vext-----[L]:</td> <td style="text-align: center;">0.078</td> <td style="text-align: center;">0.053</td> <td style="text-align: center;">67.9%</td> </tr> <tr> <td>EX Time---[s]:</td> <td style="text-align: center;">2.4</td> <td style="text-align: center;">1.7</td> <td style="text-align: center;">70.8%</td> </tr> <tr> <td>MMF-----[L/s]:</td> <td style="text-align: center;">2.763</td> <td style="text-align: center;">2.991</td> <td style="text-align: center;">108.2%</td> </tr> <tr> <td>PEF-----[L/s]:</td> <td style="text-align: center;">5.577</td> <td style="text-align: center;">5.399</td> <td style="text-align: center;">96.8%</td> </tr> <tr> <td>MEF75---[L/s]:</td> <td style="text-align: center;">4.986</td> <td style="text-align: center;">5.091</td> <td style="text-align: center;">102.1%</td> </tr> <tr> <td>MEF50---[L/s]:</td> <td style="text-align: center;">3.359</td> <td style="text-align: center;">3.536</td> <td style="text-align: center;">105.2%</td> </tr> <tr> <td>MEF25---[L/s]:</td> <td style="text-align: center;">1.437</td> <td style="text-align: center;">1.608</td> <td style="text-align: center;">111.8%</td> </tr> <tr> <td>FIVC-----[L]:</td> <td style="text-align: center;">2.494</td> <td style="text-align: center;">2.392</td> <td style="text-align: center;">95.9%</td> </tr> </tbody> </table> <div style="text-align: center; font-size: x-small;">[▲][▼]Page Up and Down</div>	Para.	FVC		2/4	Pre.	Post.	%Pre.	Vext-----[L]:	0.078	0.053	67.9%	EX Time---[s]:	2.4	1.7	70.8%	MMF-----[L/s]:	2.763	2.991	108.2%	PEF-----[L/s]:	5.577	5.399	96.8%	MEF75---[L/s]:	4.986	5.091	102.1%	MEF50---[L/s]:	3.359	3.536	105.2%	MEF25---[L/s]:	1.437	1.608	111.8%	FIVC-----[L]:	2.494	2.392	95.9%
Para.		FVC		1/4																																																																											
	Pre.	Post.	%Pre.																																																																												
FVC-----[L]:	2.671	2.493	93.3%																																																																												
FEV0.5---[L]:	1.753	1.751	99.8%																																																																												
FEV1.0---[L]:	2.276	2.247	98.7%																																																																												
FEV3.0---[L]:	2.616	2.463	94.1%																																																																												
FEV1.0/G---[%]:	85.2%	99.8%	105.7%																																																																												
FEV1.0/T---[%]:	84.2%	89.8%	106.6%																																																																												
FEV3.0/G---[%]:	97.9%	98.7%	100.8%																																																																												
FEV3.0/T---[%]:	96.8%	98.4%	101.6%																																																																												
Para.	FVC		2/4																																																																												
	Pre.	Post.	%Pre.																																																																												
Vext-----[L]:	0.078	0.053	67.9%																																																																												
EX Time---[s]:	2.4	1.7	70.8%																																																																												
MMF-----[L/s]:	2.763	2.991	108.2%																																																																												
PEF-----[L/s]:	5.577	5.399	96.8%																																																																												
MEF75---[L/s]:	4.986	5.091	102.1%																																																																												
MEF50---[L/s]:	3.359	3.536	105.2%																																																																												
MEF25---[L/s]:	1.437	1.608	111.8%																																																																												
FIVC-----[L]:	2.494	2.392	95.9%																																																																												

FVC				3/4
Para.	Pre.	Post.	%Pre.	
FIV0.5	[L]: 1.012	0.769	75.9%	
FIV1.0	[L]: 2.084	2.044	98.0%	
FIV1.0/FVC	[%]: 78.0%	81.9%	105.0%	
FIV1.0/FIVC	[%]: 83.5%	85.4%	102.2%	
PIF	[L/s]: 2.641	2.635	99.7%	
MIF50%	[L/s]: 2.344	2.011	111.3%	

[▲][▼]Page Up and Down


FVC Test Result		4/4
Diagnosis:		
Lung diseases and degree:		
Mild obstruction		
Improvement (Post):		
No improved		

[START]Retest, [▲]Page Up, [↵]Enter

Patient Test Data	ID: 000000015 SVC 2009/08/27 14:07	ID: 000000015 MVV 2009/08/27 14:07
Name: _____	Para Unit Pred. Meas. %Pred. Post. %Pre.	Para Unit Pred. Meas. %Pred. Post. %Pre.
2009/08/27 14:07	VC L 4.613 3.767 81.6% 4.015 106.5%	MVV LPM 158.484 119.640 75.4% 107.600 89.9%
ID No.: 0000000015	TV L 0.501 0.856 170.8%	RR BPM 60 49 81.6%
Gender: Male Age: 30 yrs	ERV L 1.494 1.569 105.0%	TV L 1.994 2.196 110.1%
Height: 170 cm Weight: 70 kg	IRV L 1.772 1.590 89.7%	
Race: Chinese	IC L 2.273 2.446 107.6%	
Temp.: 29 C Humi.: 40 %		
Atoms: 760 mmHg		
Test Result:		
Doctor: _____		

Post-measurement report part 1

ID: 0000000016 FVC 2009/08/27 14:27						
Para	Unit	Pred.	Meas.	%Pred.	Post.	%Pre.
FVC	L	4.613	3.671	79.5%	3.777	102.8%
FEV0.5	L	2.997	2.446	81.6%	2.387	97.5%
FEV1.0	L	3.869	3.202	82.7%	3.150	98.3%
FEV3.0	L	4.483	3.644	81.2%	3.715	101.9%
FEV1.0%G	%	83.8%	87.2%	104.0%	83.3%	95.5%
FEV1.0%T	%	0.0%	0.0%	0.0%	0.0%	-----
FEV3.0%G	%	97.0%	99.2%	102.2%	98.3%	99.0%
FEV3.0%T	%	0.0%	0.0%	0.0%	0.0%	-----
Vext	L		0.063		0.042	66.6%
EX Time	s		2.1		3.6	171.4%
MMF	L/s	4.193	3.671	87.5%	3.237	88.1%
PEF	L/s	8.400	9.094	108.2%	8.807	96.8%
MEF75	L/s	7.794	7.059	90.5%	6.631	93.9%
MEF50	L/s	4.689	4.543	96.8%	3.996	87.9%
MEF25	L/s	1.971	1.884	95.5%	1.536	81.5%
FIVC	L		3.337		3.486	104.4%
FIV0.5	L		1.514		1.141	75.3%
FIV1.0	L		3.222		2.830	87.8%
FIV1.0/FVC	%		87.7%		74.9%	85.4%
EIV1.0/FIVC	%		96.5%		81.1%	84.0%
PIF	L/s		4.943		3.843	77.7%
MIF50%	L/s		4.933		3.628	73.5%




Diagnosis: _____

Lung diseases and degree: _____

Normal

Improvement (Post): No improved

Post-measurement report part 2



Post-measurement report part 3