


AED Accessories















• www.meditech.com.cn

MEDITECH







MEDITECH INDOOR AED CABINET

Item	Image	Alarm	Details
MDA-E01			<ol style="list-style-type: none"> 1).Size:42x37x15.5cm 2).Weight:6.5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Alarm System with 1.5V battery x 3pcs
MDA-E02			<ol style="list-style-type: none"> 1).Size 42x37x15.5cm 2).Weight:6.5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).9V Alarm with remote controller
MDA-E03			<ol style="list-style-type: none"> 1).Size:40x40x20cm 2).Weight:5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Alarm System with 1.5V battery x 3pcs
MDA-E04			<ol style="list-style-type: none"> 1).Size:40x40x20cm 2).Weight:5.5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).key-activated alarm,9V battery operated 85db to 120db door alarm, louder than other units
MDA-E05			<ol style="list-style-type: none"> 1).Size:38x38x18cm 2).Weight:5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Alarm System with 1.5V battery x 3pcs
MDA-E06			<ol style="list-style-type: none"> 1).Size:40x38x20cm 2).Weight:6.5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).9V Alarm with remote controller
MDA-E07			<ol style="list-style-type: none"> 1).Size:46x46x11cm 2).Weight:7kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Optional key-activated,9V battery operated 85db to 120db door alarm, louder than other units







MDA-E08			<ol style="list-style-type: none"> 1).Size:36x36x18cm 2).Weight:5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Alarm System with 1.5V battery x 3pcs
MDA-E09			<ol style="list-style-type: none"> 1).Size 40x37x20cm 2).6.5kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Optional key-activated,9V battery operated 85db to 120db door alarm, louder than other units
MDA-E10			<ol style="list-style-type: none"> 1).Size:40x40x20cm 2).Weight:8kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Optional key-activated,9V battery operated 85db to 120db door alarm, louder than other units 6).Strobe light 7).Lock
MDA-E11			<ol style="list-style-type: none"> 1).Size:40x38x20cm 2).Weight:5.5kg 3).Plexiglas glass in front. 5).Wall-mounted 7).Strobe light with sounder
MDA-E12			<ol style="list-style-type: none"> 1).Size:40x38x20cm 2).Weight:6kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Optional key-activated,9V battery operated 85db to 120db door alarm, louder than other units 6).Strobe light 7).Lock
MDA-E13			<ol style="list-style-type: none"> 1).Size:40x40x20cm 2).Weight:7kg 3).Plexiglas glass in front. 4).Wall-mounted 5).Optional key-activated,9V battery operated 85db to 120db door alarm, louder than other units 6).Strobe light
MDA-E14			<ol style="list-style-type: none"> 1).Size:46x36x26cm 2).Weight:8kg 3).Plexiglas glass in front. 4).Alarm system 5).Wall-mounted 6).Adjustable Shelf

MEDITECH OUTDOOR AED CABINET


MDA-H1			<ol style="list-style-type: none"> 1).Size:48x38x24cm 2).Weight:10kg 3).Temperature display(220V 10A) 4).Thermostat(220v 100W) 5).Temperature rang:0-70°C 6).Alarm inside. When the door is opened the alarm will sound
MDA-H2			<ol style="list-style-type: none"> 1).Size:48x38x24cm 2).Weight:11kg 3).Temperature display(220V 10A) 4).Thermostat(220v 100W) 5).Temperature rang:0-70°C 6).Alarm inside. When the door is opened the alarm will sound
MDA-H3			<ol style="list-style-type: none"> 1).Size:46x36x26cm 2).Weight:12kg 3).Plexiglas glass in front. 4).Alarm system 5).Wall-mounted 6).Adjustable Shelf 7).Heating System
MDA-H4			<ol style="list-style-type: none"> 1).Size:48x38x24cm 2).Weight:15kg 3).Temperature display(220V 10A) 4).Thermostat(220v 100W) 5).Temperature rang:0-70°C 6).Alarm inside. When the door is opened the alarm will sound 7).Mechanic Lock 8).Light (door opened light on) (door closed light off)
MDA-H3-S			<ol style="list-style-type: none"> 1).Size:48x38x24cm 2).Weight:11kg 3).Alarm inside. When the door is opened the alarm will sound 4).Mechanic Lock <p>NOTE: this one does not have heating system</p>
MDA-H2-S			<ol style="list-style-type: none"> 1).Size 48x38x24cm 2).Weight:9.5kg 3).Alarm inside. When the door is opened the alarm will sound 4).Digital Lock <p>NOTE: this one does not have heating system</p>








MDA-H1-S			<p>1).Size:48x38x24cm 2).Weight:9kg 3).Alarm inside. When the door is opened the alarm will sound 4).Regular Lock</p> <p>NOTE: this one does not have heating system</p>
MDA-E25			<p>1).Size:40x38x20cm 2).Weight:6kg 3).Alarm inside. When the door is opened the alarm will sound 4).Lock 5).Emergency Key</p>


MEDITECH STANDING AED CABINET

MDAHD-S01			<p>1).Size 140x48x24cm 2).Material:1.5mm cold rolled steel 3).Outdoor 4).Alarm inside. When the door is opened the alarm will sound</p>
MDAHD-S02			<p>1).Size 140x45x22cm 2).Material:1.5mm cold rolled steel 3).Indoor 4).Alarm inside. When the door is opened the alarm will sound</p>
MDAHD-S03			<p>1).Size 140x48x24cm 2).Material:1.5mm cold rolled steel 3).Outdoor 4).Temperature display(220V 10A) 5).Thermostat(220v 100W) 6).Temperature rang:0-70°C 7).Alarm inside. When the door is opened the alarm will sound</p>

MEDITECH AED SIGN

Item	Image	Specification
MDA-S-1		<p>1).Size 30×20cm 2).Aluminum 3).White</p>

MDA-S-2		<p>1).Size 23×17cm 2).Aluminum 3).Green</p>
MDA-S-3		<p>1).Size 31×23cm 2).Aluminum 3).White</p>
MDA-S-4		<p>1).Size 45.7×20cm 2).Aluminum 3).Green</p>
Meditech AED BRACKET		
MDA-B-01		<p>1).Size 23x17.8x13.5cm 2).Metal(cold rolled steel) 3).Powder coating</p>
MDA-B-02		<p>1).Size 40x25x15cm 2).Metal(cold rolled steel) 3).Powder coating</p>
MDA-B-03		<p>1).Size 25x17.8x13.5cm 2).Metal(cold rolled steel) 3).Powder coating</p>
MDA-B-04		<p>1).Size 25x25x15cm 2).Metal(cold rolled steel) 3).Powder coating</p>

MDA-B-05		1).Size 19x12.3x9.5cm 2).Metal(cold rolled steel) 3).Powder coating
----------	--	---



Meditech Equipment Co., Ltd (Meditech Group)

Tel: (+86) 532 85832673 81705331

Fax: (+86) 532 81705332

Email: sales@meditech.cn

USA : (+1) 509 588 7272

Email: US@meditech.cn

Europe: (+44) 020 3239 3734

Email: UK@meditech.cn

Middle East: (+2) 02 27740444

Email: mdm@meditech.cn

Distributor :



Meditech Customer Care Team

To serve you well please kindly CC: every mail you send to meditech contact person to:

Customer Care Team: cct@meditech.cn

If you want to know the distributor of your country please contact us at :

cct@meditech.cn ,and we`ll get back to you as soon as possible.

After sale Service: service@meditech.cn